

HOUSE/HOUGH STYLE HANDLE


- Used as a connector between the suction tip and the suction tubing
- Sterile, packaged in individual peel pouch

Rhinologic

Number	Description	Qty/Box
66500992	Suction Tip Handle with control hole	10
66505009	Suction Tip Handle without control hole	10
66500990	Elongated House/Hough style handle with finger control hole, 3"	10

Otologic

Number	Description	Qty/Box
66500992	Suction Tip Handle with control	10
66505009	Suction Tip Handle without control	10
66500990	Elongated House/Hough style handle with finger control hole, 3"	10



METAL SUCTION TIP

- Applicable for ear examinations, suction of wax and other fluids, paracentesis
- The curved angle is specially designed for microscopic examination and surgery
- Quick connector suction tip handle is easy to change during surgery
- Sterile, packaged in individual peel pouch
- Available in broad variety of lumen sizes


Rhinologic

Number	Description	Color	Qty/Box
65000045	0.4 mm x 80 mm	Lilac	10
66500072	0.7 mm x 80 mm	Black	10
66500100	1.0 mm x 70 mm, 3 Fr	Brown	10
66500102	1.0 mm x 80 mm, 3 Fr	Brown	10
66500140	1.4 mm x 70 mm, 4 Fr	White	10
66500142	1.4 mm x 80 mm, 4 Fr	White	10
66500200	2.0 mm x 80 mm, 14 g, 6 Fr	Green	10
66500250	2.5 mm x 80 mm, 7.5 Fr	Orange	10
66500300	3.0 mm x 80 mm, 9 Fr	Blue	10
66500400	4.0 mm x 80 mm, 12 Fr	Yellow	10
66509100	1.0 mm x 100 mm, 3 Fr	Brown	10
66509140	1.4 mm x 100 mm, 4 Fr	White	10
66509200	2.0 mm x 100 mm, 6 Fr	Green	10
66509250	2.5 mm x 110 mm, 7 Fr	Orange	10
66509300	3.0 mm x 100 mm, 9 Fr	Blue	10
66516020	2.0 mm x 150 mm, 6 Fr	Green	10
66516030	3.0 mm x 150 mm, 9 Fr	Yellow	10
66516040	4.0 mm x 150 mm, 12 Fr	Blue	10



ROSEN-HOUSE STYLE SUCTION NEEDLE


Number	Description	Dimensions	Color	Qty/Box
66500040	Micro suction tip	0.4 mm x 70 mm, 26 G	Lilac	10
66500045	Micro suction tip	0.4 mm x 80 mm, 26 G	Lilac	10
66500070	Micro suction tip	0.7 mm x 70 mm, 22 G	Black	10
66500072	Micro suction tip	0.7 mm x 80 mm, 22 G	Black	10
66500100	Micro suction tip	1.0 mm x 70 mm, 20 G, 3 Fr	Brown	10
66500140	Micro suction tip	1.4 mm x 70 mm, 18 G, 4 Fr	White	10
66500208	Suction tip	2.0 mm x 70 mm, 14 G, 5 Fr	Green	10
66500250	Suction tip	2.5 mm x 80 mm, 7 Fr	Orange	10
66500300	Micro suction tip	3.0 mm x 80 mm, 9 Fr	Blue	10



FRAZIER/FERGUSON NASAL STYLE SUCTION

- Designed for Rhinology surgery or examination
- Suction tip made of medical grade stainless steel
- Sterile, packed in individual peel pouch
- Available with one side hole (velvet eye) at tip end to avoid biopsy effect (tissue damage)


Number	Description / Dimensions	Color	Qty/Box
66501020	2.0mm x 150mm, 6 Fr	Clear	10
66501030	3.0mm x 150mm, 9 Fr	Clear	10
66501040	4.0mm x 150mm, 12 Fr	Clear	10
66511020	2.0mm x 150mm, 6 Fr with side hole	Clear	10
66511030	3.0mm x 150mm, 9 Fr with side hole	Clear	10
66511040	4.0mm x 150mm, 12 Fr with side hole	Clear	10



V. EICKEN SINUS ANTRUM CANNULA, LONG CURVE

- FESS procedures
- Olive ball tip for atraumatic suction dissection or removal of synechia
- Pencil grip
- Sterile


Number	Description	Dimensions	Color	Qty/Box
66503020	SS Olive Tip Suction	2.0mm x 150mm, 6", long curve	Clear	10
66503025	SS Olive Tip Suction	2.5mm x 150mm, 6", long curve	Clear	10
66503030	SS Olive Tip Suction	3.0mm x 150mm, 6", long curve	Clear	10
66503035	SS Olive Tip Suction	3.5mm x 150mm, 6", long curve	Clear	10
66503040	SS Olive Tip Suction	4.0mm x 150mm, 6", long curve	Clear	10



SINGLE USE PLASTIC TIPS


- Used during Rhinology, Laryngology, and Sinus procedures
- Plastic suction tip with smoothly rounded catheter tip
- Sterile product, packed in individual peel pouch

Number	Description	Dimensions	Color	Qty/Box
66100635	Plastic Tip, curved	3.0mm x 160mm, 9 Fr	Blue	10
66100643	Plastic Tip, curved	4.0mm x 160mm, 12 Fr	White	10
66100627	Plastic Tip, curved	2.0mm x 160mm, 6 Fr	Green	10



MCGEE/SCHUKNECHT STYLE FOOT PLATE SUCTION TIP


Number	Description	Dimensions	Color	Qty/Box
66507040	Baron/Rosen step down Needle Suction	0.7 (22 G) to 0.4 mm x 77 mm tip, 26 g, shaft 22 g	Lilac	10



INVOTEC SINGLE-USE SUCTION SETS

- Sterile product, packaged in individual peel pouch

Number	Description	Qty/Box
66502070	Middle Ear Surgery Set includes: Metal Suction Tips (2 each 1.0mm, 1.4mm and 2.00mm), handle and 3.5m suction tube	1 set
66502208	Suction Kit Set includes: Metal Suction Tip (2.0mm x 70mm), handle and 2m long suction tube	1 set
66502140	Suction Kit Set includes: Metal Suction Tip (1.4mm x 70mm), handle and 2m long suction tube	1 set




MEDIBENDER

A bending instrument for Metal Suction Tips

- Stainless steel
- Non-sterile
- Autoclavable

Number	Description	Qty/Box
6066500	Medibender for use with 2 mm -1.4 mm -1 mm Ø Single sided	1
6066501	Medibender for use with 3 mm (9 FR) - 4 mm (12 FR) Ø Double sided	1



SUCTION TUBING

Number	Description	Length	Qty/Box
66153510	Suction Tubing for metal tips	PVC 3.5 mm length	5
6099996	Suction Tubing for plastic tips	PVC 3.5 mm length	5