


### BREATHE-EASY™ SPLINT

Silicone nasal septal splint with integral airway

- Provides septal support and allows nasal breathing post-operative through integral airway
- Pre-punched hole in anterior tip for easier suturing to prevent dislodgement and aid in septal stabilization
- 100% medical grade silicone for ease of insertion and removal
- Provided sterile, ready to use

Number	Description	Size	Qty/Box
20-10500	Breathe-Easy Nasal Splint	Regular	5 pair
20-10500BSP	Breathe-Easy Nasal Splint	Regular	50 pair
20-10510	Breathe-Easy Nasal Splint	X-Large	5 pair




### PRE-CUT™ SPLINT

Silicone nasal septal splint

- Provides septal support
- Pre-punched hole in anterior tip for easier suturing to prevent dislodgement and aid in septal stabilization
- 100% medical grade silicone for ease of insertion and removal
- Provided sterile, ready to use

Number	Description	Size	Qty/Box
20-10505	Pre-Cut Nasal Splint	Sterile	5 pair





### DEMPSEY SPLINT

Fluoroplastic internal nasal splint

- 100% Fluoroplastic to reduce or prevent adhesions between the septum and lateral nasal wall by stabilizing the turbinates
- Easy insertion and removal
- Provided sterile, ready to use

Number	Description	Size	Qty/Box
20-10602	Dempsey Medium	Sterile	10 pair
20-10603	Dempsey Large	Sterile	10 pair
20-10602C	Dempsey, Medium, Clear	Sterile	10 pair
20-10603C	Dempsey, Large, Clear	Sterile	10 pair




### K-SPLINT™

Designed to be used for endoscopic sinus surgery or any other procedure where middle turbinate adhesions or lateralizations are problematic

- Easy to shape and configure to patient's septum and middle turbinate anatomy
- Available in four sizes to reduce sizing time in the operating room
- Made of easy to see white malleable fluoroplastic
- Supplied sterile for immediate surgical usage
- Designed by J. Patrick Shea, M.D.

Number	Description	Size	Qty/Box
20-10640	Small	Sterile	5 pair
20-10641	Medium/Small	Sterile	5 pair
20-10642	Medium	Sterile	5 pair
20-10644	Large	Sterile	5 pair




## GUATELLA/MANTOVANI™ INTERNAL NASAL SPLINT

- Designed for functional and/or aesthetic rhinosurgery of the septal valve nose tip and nasal pyramid
- Shaped to support the tip of the nose
- Made from 100% Fluoroplastic
- Packaged sterile, ready to use

Number	Description	Thickness	Qty/Box
20-10630	Guastella/Mantovani Internal Nasal Splint	.25mm	5 pair

Reference: *Mantovani, M., Guastella, C., Mazzola, R.*  
*The Guastella/Mantovani Septal-Valve Splint: An Intranasal*  
*Biplanar and Multivectorial Orthopedic Device*  
*Plastic and Reconstructive Surgery, Vol. 106, August, 2000*



## BIVALVE SPLINTS

### Nasal Septal Splints

- 100% Fluoroplastic to reduce or prevent adhesions between the septum and lateral nasal wall
- Easy insertion and removal
- Pre-punched suture holes
- Over 30 years of clinical use
- Provided sterile and ready to use

## CUSTOM-CUT™ SPLINTS

- Choice of 100% Fluoroplastic or 100% medical grade silicone
- Allows for custom design by the surgeon
- Provided sterile and ready to use
- Five thicknesses available

Number	Size	Description	Thickness		Qty/box
20-10650	Standard	Bivalve Splint	.25mm	Sterile	5 pair
20-10650S	Standard	Bivalve Splint	.25mm	Sterile	10 ea
20-10655	Standard	Bivalve Splint	.50mm	Sterile	5 pair
20-10655S	Standard	Bivalve Splint	.50mm	Sterile	10 ea
20-10660	Large	Bivalve Splint	.25mm	Sterile	5 pair
20-10660S	Large	Bivalve Splint	.25mm	Sterile	10 ea
20-10665	Large	Bivalve Splint	.50mm	Sterile	5 pair
20-10665S	Large	Bivalve Splint	.50mm	Sterile	10 ea
20-10670	38mm x 64mm	Custom-Cut Fluoroplastic	.25mm	Sterile	5 pair
20-10675	38mm x 64mm	Custom-Cut Fluoroplastic	.50mm	Sterile	5 pair
20-10680	55mm x 75mm	Custom-Cut Silicone	1.0mm	Sterile	5 pair
20-10680S	55mm x 75mm	Custom-Cut Silicone	1.0mm	Sterile	10 ea
20-10681	55mm x 75mm	Custom-Cut Silicone	.25mm	Sterile	10 ea
20-10682	55mm x 75mm	Custom-Cut Silicone	.50mm	Sterile	5 pair
20-10682S	55mm x 75mm	Custom-Cut Silicone	.50mm	Sterile	10 ea
20-10685	38mm x 50mm	Custom-Cut Silicone	.13mm	Sterile	10 ea
20-10685B	38mm x 50mm	Custom-Cut Silicone, Blue	.13mm	Sterile	10 ea
20-10686	55mm x 75mm	Custom-Cut Silicone	1.5mm	Sterile	5 pair
20-10686S	55mm x 75mm	Custom-Cut Silicone	1.5mm	Sterile	10 ea